



TechPPE

0141 761 4700

sales@iclttech.co.uk

www.techppe.co.uk

Haika HKC02 Reusable Surgical Gown

HAIKA

Key Features

High quality finish, manufactured in the UK

Machine washable up to 73°C for 100 washes (currently in testing)

Elasticated soft touch ribbed cuffs

Thumb loops for added security

Ergonomic fit for freedom of movement

One size fits all wraparound design for enhanced coverage

Secure back neck studded adjustment

Fluid repellent PVC coating with soft comfortable feel

Stitched and high frequency top welded seams for greater protection

Easily adjusted with side tie fastening

Quick drying



Made in UK



Medical device

EN13795-1

Ordering Information

- **Product Code:** HKC00002EZ
- **Pack Quantity:** 1 single per bag
- **Case Quantity:** 10 singles per carton

EU Compliance

- Class 1 Medical Device
- EN13795-1 Standard Performance



TechPPE

0141 761 4700

sales@icltech.co.uk

www.techppe.co.uk

Haika HKC02 Reusable Surgical Gown



Storage and Maintenance

Wash Instructions

- Non-ionic detergent should be used
- Do not use bleach
- Do not iron
- Do not dry clean
- Wash at a maximum temperature of 75°C
- Hang dry to minimise creasing
- Maximum 100 washes (under testing)

Storage

- It is recommended that the garments are stored in their original packaging clear of UV sources. Following use and cleaning, garments should be stored in a cool, dry and clean environment.

Components

- **Body:** 81% PVC 19% Polyester
- **Cuff:** Polyester
- **Thumb Loop:** Braided Elastic
- **Tie:** Polyester Herringbone Weave
- **Binding:** Polyester Blend

Sizing

- **Item:** One size
- **Center Back Length (cm):** 120
- **Width (cm):** 160
- **Under Arm Length including cuff (cm):** 64.5



A design collaboration between ICL Tech Ltd and Globus Group.
Manufactured by Globus in the UK // Distributed by ICL Tech Ltd.



Supplied by ICL Tech Ltd, 1 Wellington Road, Bishopbriggs, Glasgow, G64 2SA

0141 761 4700

sales@icltech.co.uk

www.techppe.co.uk